


# Supplementary Methods

## pCFD5 cloning protocol



pCFD5 is a vector for expressing one or multiple tRNA-flanked Cas9 gRNAs under the control of the strong, ubiquitous RNA pol III promoter, U6:3.

### Prepare the backbone:

Step 1: Digest the pCFD5 plasmid with Bbs1 type-IIS restriction enzyme (NEB). We typically digest 8 µg DNA with 1 µl enzyme (10 u) and 3 µl 10x buffer in a 30 µl reaction for 2 - 4 h at 37°C.




The BbsI digest creates cohesive ends that allow seamless cloning of guide sequences.

Step 2: Run the digested plasmid on a 1% agarose gel. Cut out the backbone (9.6kb) and purify DNA using a commercial kit. Elute in 25 µl sterile dH<sub>2</sub>O and measure DNA concentration. Concentration should be 50 - 200 ng/µl, which is enough for several cloning reactions.

### Cloning one gRNA:

Cloning a single gRNA into the pCFD5 plasmid uses annealed oligos, analogous to the cloning method used for pCFD3 (Port et al. 2014). It is fast, cheap and highly efficient. Thus, many plasmids can be produced in parallel, if desired.

Step 1: Order standard desalted oligos without further purification:



The guide sequence matches the genomic target site. The length of the guide sequence is usually 20 nt, but truncations are often tolerated and can make the CRISPR system more specific (Fu et al., 2014). Note that, whereas pCFD1 - 4 require a G at the first position of the guide sequence for optimal RNA pol III transcription, there is no such constraint for pCFD5 as transcription is initiated before the 5' tRNA.

Step 2: Resuspend oligos in sterile dH<sub>2</sub>O to a concentration of 100 µM.

## Cloning one gRNA (continued):

Step 3: Set up the following phosphorylation and annealing reaction:

1 $\mu$ l	Top oligo (100 $\mu$ M stock)
1 $\mu$ l	Bottom oligo (100 $\mu$ M stock)
1 $\mu$ l	10x T4 Ligation Buffer (NEB)
6.5 $\mu$ l	ddH <sub>2</sub> O
0.5 $\mu$ l	T4 Polynucleotide Kinase (NEB)

Step 4: Incubate and anneal in a thermocycler:

37°C	30 min
95°C	5 min
ramp down to 20°C at 5°C/min	

Step 5: Set up the following ligation reaction:

X $\mu$ l	BbsI digested pCFD5 (60 ng)
1 $\mu$ l	annealed oligos diluted 1:200 in ddH <sub>2</sub> O
1.5 $\mu$ l	10x T4 Ligation Buffer (NEB)
X $\mu$ l	ddH <sub>2</sub> O
1 $\mu$ l	T4 DNA Ligase (NEB)
Total volume 15 $\mu$ l	


Step 6: Ligate 30 min at room temperature.

Step 7: Transform into competent bacteria. We use 2  $\mu$ l of the ligation reaction to transform 70  $\mu$ l chemically competent bacteria (Silver competent cells, Bioline). Let cells recover in 700  $\mu$ l 2xTY media and plate out 100  $\mu$ l on Ampicillin plates.

Step 8: Incubate plates at 37°C overnight.

Step 9: Pick two colonies and verify inserts by Sanger sequencing with primer U63seqfwd: 5'-ACGTTTATAACTTATGCCCTAAG-3'

## Cloning 2 - 4 gRNAs:



The individual guide sequences are encoded in PCR primers that are used on a pCFD5 template. The resulting PCR products are assembled with the linearised pCFD5 backbone in a single Gibson Assembly reaction.



## Cloning of 2 - 4 gRNAs (continued):

Step 2: Run the individual PCRs with the primers above and *pCFD5* as template. Use a high-fidelity polymerase according to the manufacturer's guidelines. We use the 2x Q5 hot-start kit from NEB and an annealing temperature of 61°C that increases 0.5°C per cycle until a temperature of 72°C is reached (total number of cycles: 32).

Step 3: Run the PCR products on a 1% agarose gel. The expected sizes are PCR1, 234 bp; PCR2, 204 bp; PCR3, 234 bp. Cut out the correct bands and extract the DNA using a commercial kit.

Step 4: Prepare a Gibson Assembly reaction using ~50 ng of digested backbone and an ~two-fold molar excess of each insert. The typical reaction volume is 15 µl. Incubate for 1 h at 50°C.

Step 5: Transform 2 µl of Gibson Assembly mix into 70 µl chemically competent cells using standard procedures (note that too high a proportion of Gibson Assembly mix interferes with transformation). Plate cells on Ampicillin plates and incubate at 37°C overnight.


### Recommended:

The assembly of more than two fragments by Gibson Assembly is error prone. To reduce the number of clones that need to be sequenced in these cases, we perform colony PCR to select clones with inserts that appear to be of the right size.

Step 6: Pick single colonies with a pipet tip or plastic inoculation needle and dissolve in 10 µl bacterial culture medium containing Ampicillin. Use 1 µl of this medium as a template in a diagnostic PCR using an inexpensive PCR system (e.g. NEB Quick load Taq 2x master mix) and the following primers:

U63seqfwd: 5'-ACGTTTTATAACTTATGCCCCCTAAG-3'

pCFDseqrev: 5'-GCACAATTGTCTAGAATGCATAC-3'



An example of a colony PCR from a Gibson Assembly of four gRNAs and *pCFD5*. Colonies marked by an asterisk had the approximate size predicted for the desired product and were analysed by sequencing.


Step 7: Examine PCR products on a 1% agarose gel. Colonies that give a product of the correct size (three gRNAs, 1013 bp; four gRNAs, 1186 bp) should be cultivated in 1 - 5 ml media overnight and miniprep DNA analysed by Sanger sequencing with primer U63seqfwd: 5'-ACGTTTTATAACTTATGCCCCCTAAG-3' and pCFDseqrev 5'-GCACAATTGTCTAGAATGCATAC-3'.

## Cloning of 5 - 6 gRNAs:

Cloning of more than four gRNAs follows the same logic as cloning three or four gRNAs. However, in our hands directly assembling more than four fragments with Gibson Assembly is inefficient. We therefore fuse individual PCR products to reduce the number of inserts to no more than three before using them in a Gibson Assembly reaction.

Cloning 5 - 6 gRNAs (continued):


I. Generate PCR products encoding all guide sequences



II. Fuse adjacent fragments by extension PCR



III. Assemble up to three fragments with the linear backbone to produce the final plasmid



Follow the protocol above for 2 - 4 gRNAs using additional internal primer pairs to generate 4 fragments (for 5 gRNAs) or 5 fragments (for six gRNAs).

Once the individual PCR products have been purified (Step 3), assemble two short adjacent fragments into one longer fragment. Use 300 ng of each PCR product in a PCR reaction mix without primers that has 15 cycles and an annealing temperature of 61°C. The overlapping ends of the two fragments will prime extension of each fragment. Add distal primers (e.g., PCR1fwd and PCR2rev) and run PCR for an additional 20 cycles with an annealing temperature of 61°C. Purify the desired product from an agarose gel. Once you have reduced the number of inserts to three, follow the protocol above for Gibson Assembly and identification of correct inserts. In Step 7, the expected product sizes for five and six gRNAs are 1360 bp and 1534 bp, respectively.

References:

Fu Y, Sander JD, Reyon D, Cascio VM, Joung JK. Improving CRISPR-Cas nuclease specificity using truncated guide RNAs. *Nat Biotechnol.* 2014 Mar;32(3):279-84. doi: 10.1038/nbt.2808. Epub 2014 Jan 26.

Port F, Chen HM, Lee T, Bullock SL. Optimized CRISPR/Cas tools for efficient germline and somatic genome engineering in *Drosophila*. *Proc Natl Acad Sci U S A.* 2014 Jul 22;111(29):E2967-76. doi: 10.1073/pnas.1405500111. Epub 2014 Jul 7.