



## Cloning of one gRNA (continued):

Step 3: Set up the following phosphorylation and annealing reaction:

1 $\mu$ l	Top oligo (100 $\mu$ M stock)
1 $\mu$ l	Bottom oligo (100 $\mu$ M stock)
1 $\mu$ l	10x T4 Ligation Buffer (NEB)
6.5 $\mu$ l	ddH <sub>2</sub> O
0.5 $\mu$ l	T4 Polynucleotide Kinase (NEB)

Step 4: Incubate and anneal in a thermocycler:

37°C	30 min
95°C	5 min
ramp down to 20°C at 5°C/min	

Step 5: Set up the following ligation reaction:

X $\mu$ l	Bbs1 digested pCFD6 (use 60 ng)
1 $\mu$ l	annealed oligos diluted 1:200 in ddH <sub>2</sub> O
1.5 $\mu$ l	10x T4 Ligation Buffer (NEB)
X $\mu$ l	ddH <sub>2</sub> O
1 $\mu$ l	T4 DNA Ligase (NEB)
Total volume 15 $\mu$ l	


Step 6: Ligate 30 min at room temperature.

Step 7: Transform into competent bacteria. We use 2  $\mu$ l ligation to transform 70  $\mu$ l chemically competent bacteria (Silver competent cells, Bionline), let cells recover in 700  $\mu$ l 2xTY media and plate out 100  $\mu$ l.

Step 8: Plate on Ampicillin plates and incubate at 37°C overnight.

Step 9: Pick two colonies and verify inserts by Sanger sequencing with primer pCFD6seqfwd: 5'-GTAGACATCAAGCATCGGTGG-3'


## Cloning of 2 - 4 gRNAs:



The individual guide sequences are encoded in PCR primers that are used on a pCFD6 template. The resulting PCR products are assembled with the linearised pCFD6 backbone in a single Gibson Assembly reaction.

## Cloning of 2 - 4 gRNAs (continued):

Below is illustrated how to clone four gRNAs from three overlapping PCR products. Two or one PCR products can be used to clone, respectively, three or two tRNA-flanked gRNAs.



Step 1: Order the oligos illustrated below:

PCR1fwd: guide sequence 1  
CGGCCCGGGTTCGATTCCCGGCCGATGCANNNNNNNNNNNNNNNNNNNNNNGTTTCAGAGCTATGCTGGAAAC

PCR1rev: guide sequence 2 (rev comp.)  
NNNNNNNNNNNNNNNNNNNNNTGCACCAGCCGGGAATCGAACC

PCR2fwd: guide sequence 2  
NNNNNNNNNNNNNNNNNNNNNGTTTCAGAGCTATGCTGGAAAC

PCR2rev: guide sequence 3 (rev comp.)  
NNNNNNNNNNNNNNNNNNNNNTGCACCAGCCGGGAATCGAACC

PCR3fwd: guide sequence 3  
NNNNNNNNNNNNNNNNNNNNNGTTTCAGAGCTATGCTGGAAAC

PCR3rev: guide sequence 4 (rev comp.)  
ATTTTAACCTTGCTATTTCTAGCTCTAAAACNNNNNNNNNNNNNNNNNNNNNTGCACCAGCCGGGAATCGAACC

Step 1: We use standard desalted oligos. Although oligos of this length can contain some mistakes without further purification, we feel that the cheaper price and faster turn-around time compensates for the additional clones that have to be sequenced. As described above, the first position of the guide sequence can be any base and the length of the guide sequence is usually 20 nt. However, truncations to the guide sequence are often tolerated and can increase specificity of mutagenesis.

Step 2: Run the individual PCRs with the primers above and pCFD6 as template. Use a high-fidelity polymerase according to the manufacturer's guidelines. We use the 2x Q5 hot-start kit from NEB and an annealing temperature of 61°C that increases 0.5°C per cycle until 72°C (total number of cycles: 32).

Step 3: Run the PCR products on an 1% agarose gel. The expected sizes are PCR1, 233 bp; PCR2, 204 bp; PCR3, 234 bp. Cut out the correct bands and extract the DNA using a commercial kit.

## Cloning of 2 - 4 gRNAs (continued):

Step 4: Prepare a Gibson Assembly reaction using ~50 ng of digested backbone and an ~two-fold molar excess of each insert. The typical reaction volume is 15 µl. Incubate for 1 h at 50°C.

Step 5: Transform 2 µl of Gibson Assembly mix into 70 µl chemically competent cells using standard procedures (note that too high a proportion of Gibson Assembly mix interferes with transformation of competent cells). Plate cells on Ampicillin plates and incubate at 37°C overnight.

### Recommended:


The assembly of > two fragments by Gibson Assembly is error prone. To reduce the number of clones that need to be sequenced in these cases, we perform colony PCR to select clones with inserts of the right size.

Step 6: Pick single colonies with a pipet tip or plastic inoculation needle and dissolve in 10 µl bacterial culture medium containing Ampicillin. Use 1 µl of this medium as a template in a diagnostic PCR using an inexpensive PCR system (e.g. NEB Quick load Taq 2x master mix) and the following primers:

fwd: GTAGACATCAAGCATCGGTGG

rev: TTAGAGCTTTAAATCTCTGTAGGTAG

Step 7: Examine PCR products on an 1% agarose gel. Colonies that give product of the correct size (three gRNAs, 1013 bp; four gRNAs, 1186 bp) should be cultivated in 1-5 ml media overnight and miniprep DNA analysed by Sanger sequencing with primer pCFD6seqfwd: 5'-GTAGACATCAAGCATCGGTGG-3'.



An example of a colony PCR from a Gibson Assembly of four gRNAs and pCFD5 (pCFD6 gives similar results). Colonies marked by an asterisk had the approximate size predicted for the desired product and were analysed by sequencing.


## Cloning of 5 - 6 gRNAs:

Cloning of more than four gRNAs follows the same logic as cloning three or four. However, in our hands directly assembling more than four fragments with Gibson Assembly is inefficient. We therefore fuse individual PCR products to reduce the number of inserts to no more than three before using them in a Gibson Assembly reaction.

continued on next page...

Cloning of 5 - 6 gRNAs (continued):


I. Generate PCR products encoding all guide sequences



II. Fuse adjacent fragments by extension PCR



III. Assemble up to three fragments with the linear backbone to the final plasmid



Follow the protocol above for 2-4 gRNAs using additional internal primer pairs to generate 4 fragments (for 5 gRNAs) or 5 fragments (for six gRNAs).

Once you have purified the individual PCR products (Step 3), assemble two short adjacent fragments into one longer fragment. Use 300 ng of each PCR product in a PCR without primers that has 15 cycles with an annealing temperature of 61°C. The overlapping ends of the two fragments will prime extension of each fragment. Add distal primers (e.g. PCR1fwd and PCR2rev) and run for an additional 20 cycles with an annealing temperature of 61°C. Once you have reduced the number of inserts to three, follow the protocol above for Gibson Assembly and identification of correct inserts. In Step 7, the expected product sizes for five and six gRNAs are 1360 bp and 1534 bp, respectively.

Noticed a mistake in this protocol? Please contact Phillip at [f.port@dkfz.de](mailto:f.port@dkfz.de)